

# Drug Abuse Test

are simple, one step test for the detection of drugs abuse in urine

Di Distribusikan Oleh :  
PT. MEDISTA UTAMA, Jakarta

Abbreviation	Drug Name	Cutoff Levels
AMP	Amphetamines	1000 ng/ml
ACE	Acetaminophen	500ng/ml
BAR	Barbiturates	300 ng/ml
BZO	Benzodiazepines	300 ng/ml
BUP	Buprenorphine	10ng/ml
COC	Cocaine	300 ng/ml
COT	Cotinine	200ng/ml
mAMP/MET	Methamphetamine	1000 ng/ml
MQL	Methaqualone	300ng/ml
EDDP	EDDP	100ng/ml
FYL	Fentanyl	200ng/ml
KET	Ketamine	1000ng/ml

feature both single and multiple panel test providing additional flexibility for targeted testing or the simultaneous detection of up to 12 drugs with a single urine sample

## Features and Benefits :

- \* Just one step
- \* Inbuilt procedural control
- \* Rapid
- \* Room temperature storage
- \* Multiple drugs measure



Urine

2°C-30°C



25 Tests/Kit, 20 Tests/kit, 1Test/kit



iSun Bio - Tech




Manufacture by:

JOYSBIO (Tianjin) Biotechnology Co., Ltd  
No. 220, Teda Dongying Rd., Tianjing City China 300457  
Telp : 86-02265372699 Fax : 86-022-65378699  
Website : www.joysbio.com

## iSun Drug Abuse Test Device (AMP/THC/MOP/MET)

Abbreviation	Drug Name	Cutoff Levels
K2	Spice (Synthetic Marijuana)	50ng/ml
MDMA	Ecstasy	500 ng/ml
MOP	Morphine	300 ng/ml
MTD	Methadone	300 ng/ml
MET	Methamphetamine	1000 ng/ml
OPI	Opiates	2000 ng/ml
OXY	Oxycodone	100 ng/ml
PCP	Phencyclidine	25 ng/ml
PPX	Propoxyphene	300 ng/ml
TCA	Tricyclic Antidepressant	1000 ng/ml
THC	Tertahydrocannabinol	50 ng/ml
TML	Tramadol	100ng/ml



Minimal training required  
Confirmation of correctly performed test  
Result in 5 minutes  
Ready to use, no need to refrigerate  
Simultaneous detection of several drugs